



## CellaVision™ DM96

*Automated Digital Cell Morphology*



# Efficiency, Proficiency, Connectivity and Collaboration

Automatic cell location and pre-classification improves resource utilization, the quality of results and employee satisfaction. One particular labor benefit is that highly skilled medical technologists are able to spend more time on difficult cases that require careful analysis and assessment. The ability to archive images enables hospitals to look at previous cell images for long term case management.

## Efficiency

Automatic cell location and pre-classification, along with unique cell views, reduce the time spent performing differentials, training new techs and monitoring proficiency.

## Proficiency


CellaVision DM96 allows different technologists to perform differentials at the same high level of consistency. The system provides traceability of results at the individual cell level.

## Connectivity

Multiple Remote Review Stations enable supervisory review, inter/intra-lab comparison and the opportunity for real-time collaboration with Pathology.

## Collaboration

As a result of network integration, the CellaVision DM96 will provide new value to clinicians, saving them time with remote, real-time access to a patient's image history.



*Images from the pre-classified cells can be magnified for detailed analysis.*

The DM96 enables laboratories to analyze blood and body fluids, including cerebrospinal, synovial and pleural fluids, in a safe, efficient and comfortable way.

Break the cycle—find the perfect workflow. The DM96 is a digital morphology system for performing:


- Automatic identification of WBCs in peripheral blood and body fluids smears
- Pre-characterization of RBC morphology
- Platelet estimation

The system's automation significantly reduces differential turnaround time, while increasing collaboration between consultants.



# Standardized Results

Utilizing CLSI guidelines in conjunction with a neural network trained on thousands of cells identified by a group of global experts, the CellaVision DM96 offers an opportunity to improve efficiency and consistency. With continued cost constraints and the shortage of skilled technologists, especially in differential analysis, the CellaVision DM96 permits technologists to use their skills to the best benefit of the laboratory.



## Laboratory Benefits

Maintain high quality by providing standardization and quick access to differential images

Access an on-board, customized, reference cell library and view cell classes side by side

Complete traceability of results at the individual cell level

Permanently store images - access a patient's image history

Remotely share images via tele-hematology. The CellaVision DM96 can share databases with other DM analyzers, installed locally or at a remote location, allowing centralized database management. By using the CellaVision Remote Review Software it is possible to transfer digital images and results within and between laboratories. The software allows for strengthened competence, qualified review and shorter turnaround times for complicated patient cases.

## Find results quickly and easily using multiple criteria

|           | Order ID | Patient ID | First Name | Last Name  | Analyzed         | Signed by |
|-----------|----------|------------|------------|------------|------------------|-----------|
| + (US)    | 20       | 1225768925 | John       | Doe        | 2004-02-02 12:08 |           |
| (US)      | 19       | 123476598  | Michael    | Hall       | 2004-02-02 12:06 |           |
| (US) (US) | 18       | 981234567  | Josie      | Clearmount | 2004-02-02 12:04 |           |
| (US)      | 16       | 987612345  | William    | Smith      | 2004-02-02 12:01 |           |
| (US) (US) | 15       | 987651234  | Claire     | Swanson    | 2004-02-02 11:59 |           |
| + (US)    | 13       | 987654312  | Joe        | Ramone     | 2004-02-02 11:56 |           |
| (US) (US) | 12       | 987654321  | John       | Smith      | 2004-02-02 11:54 |           |
| (US)      | 11       | 123456789  | Tom        | Johnson    | 2004-02-02 11:52 |           |


Ability to search the database by 15 criteria including: doctor, technologist and comments

# Diff at a Glance

Quickly screen for abnormalities, allowing fast confirmation of the CBC analyzer's results. When classifications need to be changed, cells can be edited into the appropriate cell class.



Enhance blood film review and make more efficient use of experienced morphologists' time. Screen for abnormalities to allow fast confirmation of the CBC analyzer's results.



Confirms cell counter results in seconds





## WBC Cells

- Cells are grouped with like cells to facilitate verification
- Cells are pre-classified into 18 classes:
  - **Leukocytes:** Segmented Neutrophils, Band Neutrophils, Eosinophils, Basophils, Lymphocytes, Monocytes, Blasts, Promyelocytes, Myelocytes, Metamyelocytes, Variant Lymphocytes, Plasma Cells, Unidentified
  - **Non-Leukocytes:** Smudge Cells, Artifacts, Giant Platelets, NRBCs, Platelet Clumps
- Handling leukopenic samples is more efficient due to the ability to merge cells from multiple slides on a single patient

| WBC                        | Count      | %          | x10e9/L    |
|----------------------------|------------|------------|------------|
| • Unidentified             | -          | -          | -          |
| • Band neutrophil          | 3          | 2.9        | 0.1        |
| • Segmented neutrophil     | 62         | 59.6       | 1.2        |
| • Eosinophil               | -          | -          | -          |
| • Basophil                 | -          | -          | -          |
| • Lymphocyte               | 20         | 19.2       | 0.4        |
| • Monocyte                 | 16         | 15.4       | 0.3        |
| • Promyelocyte             | -          | -          | -          |
| • Myelocyte                | -          | -          | -          |
| • Metamyelocyte            | -          | -          | -          |
| Immature eosinophil        | -          | -          | -          |
| Immature basophil          | -          | -          | -          |
| Promonocyte                | -          | -          | -          |
| Prolymphocyte              | -          | -          | -          |
| • Blast (no lineage spec)  | -          | -          | -          |
| • Lymphocyte, variant form | 3          | 2.9        | 0.1        |
| • Plasma cell              | -          | -          | -          |
| Large granular lymphocyte  | -          | -          | -          |
| Hairy cell                 | -          | -          | -          |
| Sézary cell                | -          | -          | -          |
| Other                      | -          | -          | -          |
| <b>Total</b>               | <b>104</b> | <b>100</b> | <b>2.0</b> |

## RBC Cells

- RBC Morphology is automatically pre-characterized for polychromasia, hypochromasia, aniso-, micro-, macrocytosis and poikilocytosis
- Multiple RBC fields can be scanned as if on a virtual microscope
- A dynamic micrometer facilitates additional RBC measurements
- Provides functionality for performing platelet estimates



| RBC             | 1 | 2 | 3 | 4 | %   |
|-----------------|---|---|---|---|-----|
| • Polychromasia | 1 | 1 | 1 | 1 | 0.6 |
| • Hypochromasia | 1 | 1 | 1 | 1 | 0.0 |
| • Anisocytosis  | 1 | 1 | 1 | 1 | 4.9 |
| • Microcytosis  | 1 | 1 | 1 | 1 | 4.9 |
| • Macrocytosis  | 1 | 1 | 1 | 1 | 1.8 |
| • Platyrrhysis  | 1 | 1 | 1 | 1 | 3.7 |
| Target cells    | 1 | 1 | 1 | 1 | 0.0 |
| Spherocytes     | 1 | 1 | 1 | 1 | 0.0 |
| Helmet cells    | 1 | 1 | 1 | 1 | 0.0 |
| Sickle cells    | 1 | 1 | 1 | 1 | 0.0 |
| Schistocytes    | 1 | 1 | 1 | 1 | 0.0 |

Comment

21. Auer Rods seen

Standard comments

| Code | Comment                    |
|------|----------------------------|
| 20   | Dohle bodies Present       |
| 21   | Auer Rods seen             |
| 22   | Moderate hypersegmentation |

Comment types

General

WBC

RBC

PLT

OK Cancel

## Easy Access to Powerful Features

- Add pre-coded or free text comments to any slide, cell class or specific cell
- Attach any number of cells and send them via email
- Customize the coloration and brightness of the cell images in your personal profile

# Remote review networking for clinicians - utilizing the power of integration

**Enables sharing of images for consultation, any time anywhere.**

The DM96 provides your hematology laboratory with the tools to serve the needs of multiple facilities, while providing added value to the various departments within your organization.



**Pathology**



**Oncology**



**CellaVision DM96 - The heart of your network sending information where and when it is needed**



**Pediatrics**

**Utilizes your existing LAN, WAN and VPN infrastructure**

## **Value-added:**

**For Clinicians** - Reduce time for consultation. Clinicians and pathologists save time by accessing their patient's images remotely. Morphology expertise can be centralized, resulting in staffing flexibility at remote locations.

**For Technologists** - Save time searching for requested smears. Multiple systems can be run from a central location eliminating sample transportation and reducing TAT. In real time, multiple cell images can be easily e-mailed to colleagues anywhere in the world without interrupting your workflow.


# Optional WAM System – enables results-processing / reflexing rules and seamless data integration



**Sysmex WAM™ (Work Area Manager)** - Designed for the high volume hematology laboratory, WAM delivers an advanced specimen process management system.

- Consistent and standardized rules
- Consolidates complex rules related to the criticality of results
- Patient delta-checking – two years of on-board data

A specimen-centric approach empowers WAM to deliver seamless integration with your Sysmex Intelligent Automation solution. WAM's rules-based, expert technology gives you unprecedented efficiency and control over workflow for advanced process control.



Sysmex WAM connects your remote labs using one central repository for hematology results. This gives your lab both control of the information and distributed processing capability. WAM easily interfaces to your LIS for demographics, orders results, repeats or add-on tests. Even if the LIS is down, WAM keeps your hematology department up and running.



# Optional Functionality

Now body fluids and blood smears can be analyzed on the CellaVision DM96 at the same time, either batched in blue or orange magazines or in a random access mode when interfaced to your Sysmex WAM or your LIS.

## Body Fluid Option

### Key features & benefits

- Designed to utilize standard cyto-centrifuged slides
- Pre-classification of seven WBC classes
- View cell classes side by side or all cells in a full screen view
- Digital scan of entire sample area
  - Navigation in a digital sample overview
  - Available in 10x or 50x magnification
- Regions of Interest:
  - Tag interesting areas for Pathology review and for educational purposes
  - Export your tagged areas into presentations and educational material
- Add pre-coded or free text comments to any slide, cell class or specific cell
- Permanently store images — access a patient's image history




*Images from the pre-classified cells can be magnified for detailed analysis.*

## Digital Scanning Option\*

### Key features & benefits

- Enhance, annotate and share digital slides with your students or colleagues
- Digital scan of desired sample area:
  - Navigation in a digital sample overview
  - Available in 10x or 50x magnification
- Regions of Interest (ROI):
  - Tag interesting areas for review or for educational purposes
  - Export your ROI into presentations and educational material
  - Long term storage of all images




Digitize and annotate desired areas of interesting specimens within hematology, pathology and cytology. Enable your laboratory to utilize the power of telepathology.



# Enhance Functionality

Enables the CBC results stored in WAM or LIS to be populated onto the main screen  
(Peripheral Blood Samples Only)



## Archiving Settings by Database

- Each database established may have its own specific archive settings
- Creating a new database now allocates optimal 20 GB
- Database query updated to address:
  - Body Fluids
  - Extracts the serial number of the analyzer on which the slide was scanned
- Enabled for SNCS™ monitoring

DM System Software v. 2.1 is included with the purchase of the Body Fluid Application

# Application Specification

## Peripheral Blood



|                                     |  |
|-------------------------------------|--|
| Slide Methods (Wedged)              | Automated slide maker and stainers<br>HemaPrep® /MiniPrep® automated<br>bloodsmearing device<br>Manual smears  |
| Stains                              | Romanowsky stains<br>(May Grünwald/Wright/Giemsa)  |
| Number of Cells Counted             | User Definable   |
| Quality Control                     | Built-in QC module for verification of<br>the cell location accuracy   |
| Throughput - Peripheral Blood       | Up to 35 slides/hour for complete<br>differential (100 WBC+RBC+PLT)<br>Up to 60 slides/hour (RBC and/or<br>PLT only)   |
| Digital Slides                      | Up to 30 slides/hour for 10x10 mm in 10x<br>Up to 2 slides/hour for 10x10 mm in<br>10x + 50x   |
| Slide Image Size - Peripheral Blood | 100 WBCs: ~5 MB<br>100 WBCs + RBC: ~6 MB<br>100 WBCs + RBC + PLT: ~6 MB  |
| Digital Slides                      | 10x10 mm in 10x: ~45 MB<br>10x10 mm in 10x + 50x: ~500 MB  |
| Result Parameters                   | WBC pre-classification: Segmented and<br>Band Neutrophils, Eosinophils,<br>Basophils, Lymphocytes, Monocytes,<br>Blast Cells, Promyelo cytes,<br>Myelocytes, Metamyelocytes, Variant<br>Lymphocytes, Plasma Cells and<br>Unidentified<br>Non WBC pre-classification: Smudge,<br>Artifacts, Giant Platelets, Platelet<br>Clumps, Erythroblasts (NRBC)<br>RBC pre-characterization: Automated<br>pre-characterization of aniso-, micro-<br>and macrocytosis, polychromasia,<br>hypochromasia and poikilocytosis is<br>performed in an overview image<br>corresponding to eight high power<br>fields (100x)<br>PLT estimate: The graphical user<br>interface allows manual estimation of<br>the PLT concentration, based on eight<br>high power fields (100x) |



# Application Specification

## Body Fluids

|                               |   |
|-------------------------------|---|
| Slide Preparation Methods     | Standard cytocentrifuge preparation<br>Default settings for<br>Shandon™/Wescor™/Statspin™   |
| Stains                        | Romanowsky stains<br>(May Grünwald/Wright/Giemsa)   |
| Number of Cells Counted       | User definable  |
| Quality Control               | Built-in QC module for verification of<br>the cell location accuracy  |
| Throughput                    | Based on 6 mm sample area<br>Up to 25 slides/hour for differential<br>(100 WBCs + 10x)<br>Up to 7 slides/hour for differential<br>(100 WBCs + 10x + 50x)  |
| Slide Image Size              | Based on 6 mm sample area<br>100 WBCs: ~5 MB<br>100 WBCs + 10x: ~10 MB<br>100 WBCs + 10x + 50x: ~150 MB   |
| Result Parameters             | WBC pre-classification: Neutrophils,<br>Eosinophils, Lymphocytes,<br>Macrophages (including Monocytes),<br>Other (Basophils, Lymphoma cells,<br>Atypical lymphocytes, Blasts and<br>Tumor cells) and Unidentified.<br>Non WBC pre-classification: Smudge<br>cells and Artefacts |
| Optional Software             | CellaVision® Body Fluid Application<br>CellaVision® Remote Review Software  |
| Recommended PC Specifications | 64 MB graphics RAM with Open<br>GL 1.2 support<br>Ethernet adapter 10/100 Mbps<br>512 MB RAM<br>100 MB free disk space<br>CPU Pentium® IV CellaVision<br>Competency Software  |
| Accessories                   | Barcode labeled slide magazines<br>Immersion oil<br>QC barcode labels   |





# Instrument Specifications

|   |  |
|---|--|
| System Components                               | PC with Windows XP<br>Slide Scanning Unit<br>CellaVision® DM Software<br>(U.S. Patent No. 6268611 and 6341180)   |
| Electrical Specifications All System Components | Voltage input: (115 VAC for US)<br>Current input: (8A for US)  |
| Size (W x D x H)                                | 530 x 600 x 630 mm<br>20.9 x 23.6 x 24.8 inches  |
| Slide Handling                                  | Continuous feed, walk away; up to<br>96 slides (8 magazines x 12 slides<br>/ magazine)<br>Requires barcode labeled slides with<br>clipped/round corners  |
| Oil Dispensing                                  | Automatic  |
| Archiving of Results and Images                 | Supported media: CD-R/CD-RW<br>and LAN   |
| Storage Capacity                                | Primary storage: On local hard drive<br>20 GB<br>Secondary storage: Unlimited when<br>transferred to external storage<br>media via LAN   |
| Printer Support                                 | Laser/inkjet printers supported by<br>Windows XP Communications<br>Bi-directional LIS support, ASTM<br>Ethernet 10/100 Mbps<br>E-mail  |
| Archiving Settings by Database                  | Each Database established may have<br>it's own specific archive settings<br>Creating a new database now<br>allocates optimal 20 GB<br>DB Query updated to address:<br>- Body Fluids<br>- Extracts the serial number of<br>the analyzer on which the slide<br>was scanned<br>Enabled for SNCS™ monitoring |
| Accessories                                     | Barcode labeled slide magazines<br>Immersion oil<br>QC barcode labels  |



**SYSMEX CORPORATION**  
1-5-1 Wakinohama-Kaigandori, Chuo-ku,  
Kobe 651-0073, Japan  
Phone 81-78-265-0521 • Fax 81-78-265-0530  
[www.sysmex.co.jp](http://www.sysmex.co.jp)

**SYSMEX AMERICA, INC.**  
One Nelson C. White Pkwy, Mundelein, IL 60060, U.S.A.  
Phone 1-847-996-4500 • Fax 1-847-996-4397  
[www.sysmex.com/usa](http://www.sysmex.com/usa)

**SYSMEX CANADA, INC.**  
5045 Oribitor Drive, Building 9, Suite 401  
Mississauga, ON L4W 4Y4  
Tel. 1-905-366-7900 Fax 1-905-366-7899  
[www.sysmex.ca](http://www.sysmex.ca)

**SYSMEX DO BRASIL INDÚSTRIA E COMÉRCIO LTDA**  
Rua Professor Algaycyr Munhoz Mader  
3775 - Bloco C - CIC Curitiba - Parana - Brasil - CEP 81350-010  
Phone 55-41-246-1314 • Fax 55-41-248-3180  
[www.sysmex.com/usa](http://www.sysmex.com/usa)

**SYSMEX EUROPE GMBH**  
Bornbarch 1, 22848 Norderstedt, Germany  
Phone 49-40-52726-0 • Fax 49-40-52726-100  
[www.sysmex-europe.com](http://www.sysmex-europe.com)

**SYSMEX DEUTSCHLAND GMBH**  
Bornbarch 1, 22848 Norderstedt, Germany  
Phone 49-40-5341020 • Fax 49-40-5232302  
[www.sysmex.de](http://www.sysmex.de)

**SYSMEX UK LTD.**  
Sysmex House, Garamond Drive, Wymbush  
Milton Keynes, Buckinghamshire, Mk8 8DF, U.K.  
Phone 44-870-9029270 • Fax 44-870-9029211  
[www.sysmex.co.uk](http://www.sysmex.co.uk)

**SYSMEX SINGAPORE PTE LTD.**  
2 Woodlands Sector 1, #01-06 Woodlands Spectrum  
Singapore 738068, Singapore  
Phone 65-6221-3629 • Fax 65-6221-3687  
[www.sysmex.com.sg](http://www.sysmex.com.sg)

**SYSMEX SHANGHAI LTD.**  
5th Floor, New Shanghai International Tower  
360 Pudong South Road  
Shanghai, 200120 China  
Phone 21-6886-3300 • Fax 21-6886-3355  
[www.sysmex.com.cn](http://www.sysmex.com.cn)

