



VES-MATIC Cube 200

*Automated Erythrocyte
Sediment Rate (ESR) Analyzer
ESR made simple, from Primary EDTA tubes*



Product Highlights

VES-MATIC Cube 200 Specifications

Accurate and Dependable

- Results normalized to a 60-minute Westergren method
- Standardized results
 - Consistent across all technologists and all shifts
- Minimal sample volume (1.5 mL)
 - Limits the need for second draws
- Traceability of all tested samples into dedicated racks
 - Easy, documented location of tested samples

Operational Efficiency

- Hematology sample racks are placed directly on the system
 - No technologist set-up time required
- Primary EDTA tube testing
 - No pipetting into ESR tubes
- Multi-wave length optical detection
 - No additional consumables or reagents to purchase
- Complete walk-away operation
 - Eliminates technologist's need to time the test or calculate sediment rate

Optical Detection Technology

The blood obtained in the lavender top (EDTA) test tube, typically used for hematology testing, is accurately mixed by the instrument. The samples then remain at rest for a predetermined amount of time to allow sedimentation to occur.

Through analog sensors, the instrument automatically determines the sedimentation level of erythrocytes. Subsequently, the sed rate is calculated and then automatically printed or shown on the display.

Throughput

- Up to 180 samples per hour
- Turnaround Time (TAT)
 - 26 minutes to first test result
 - Subsequent results every 18 seconds

Sampler

- 50 samples on board processing

Sample Volume

- Minimum sample volume: 1.5 mL
- Maximum sample volume: 4.0 mL

Communication Support

- Bi-directional LIS support, ASTM, via RS232C or TCP/IP conversion
- Tablet PC support, Compact Flash

Data Storage

- 10,000 samples

Printer

- Alpha-numeric with thermal paper (2 ¼") prints 36 characters per line.

Dimensions / Weight w x d x h [inches]

- Main Unit:
25 ½" x 25 ¾" x 26 ¼" / 202 pounds

Sysmex Corporation
1-5-1 Wakinohama-Kaigandori,
Chuo-ku, Kobe 651-0073, Japan
Tel. +81 (78) 265-0521
Fax +81 (78) 265-0530
www.sysmex.co.jp

Sysmex America, Inc.
One Nelson C. White Pkwy,
Mundelein, IL 60060, U.S.A.
Tel. +1 (847) 996-4500
Fax +1 (847) 996-4397
www.sysmex.com/us

Sysmex Canada, Inc.
5045 Orbitor Drive,
Building 9, Suite 401,
Mississauga, ON L4W 4Y4, Canada
Tel. +1 (905) 366-7900
Fax +1 (905) 366-7899
www.sysmex.ca

Sysmex do Brasil Indústria e Comércio Ltda
Rua Joaquim Nabuco, 615 - Bairro Cidade Jardim,
São José dos Pinhais
Paraná - Brasil - CEP 83040-210
Tel. +55 (41) 2104-1314
Fax +55 (41) 2104-1300
www.sysmex.com.br

