



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2024-08-11	Quality Control Data Due Date #1:	01-Jul-24	Quality Control Data Due Date #2:	11-Aug-24
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	41411301	41411302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.074	0.043 - 0.105	0.312	0.225 - 0.400
RBC-BF (10 ⁶ /μL)	0.025	0.020 - 0.029	0.077	0.067 - 0.088
TC-BF# (10 ³ /μL)	0.074	0.043 - 0.105	0.312	0.225 - 0.400
MN# (10 ³ /μL)	0.019	0.006 - 0.032	0.082	0.045 - 0.119
PMN# (10 ³ /μL)	0.055	0.025 - 0.084	0.230	0.157 - 0.304
MN%	26.1	10.2 - 42.0	26.3	15.5 - 37.0
PMN%	73.9	58.4 - 89.4	73.7	63.4 - 84.1
WBC ¹ (10 ³ /μL)	0.074	0.043 - 0.105	0.312	0.225 - 0.400

Please refer to the package insert for the handling of XN CHECK BF.

¹WBC for analysis in Blood Bank mode.

Distributed by:
Sysmex America, Inc.
Lincolnshire, IL

350602
Rev. 3 25-Jan-23
May-24
1



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2024-08-11	Quality Control Data Due Date #1:	01-Jul-24	Quality Control Data Due Date #2:	11-Aug-24
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	41411301	41411302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.076	0.044 - 0.108	0.314	0.226 - 0.402
RBC-BF (10 ⁶ /μL)	0.025	0.020 - 0.030	0.078	0.067 - 0.089
TC-BF# (10 ³ /μL)	0.076	0.044 - 0.108	0.314	0.226 - 0.402
MN# (10 ³ /μL)	0.020	0.009 - 0.032	0.084	0.046 - 0.122
PMN# (10 ³ /μL)	0.056	0.026 - 0.086	0.230	0.156 - 0.303
MN%	26.6	13.9 - 39.4	26.9	15.8 - 37.9
PMN%	73.4	58.0 - 88.8	73.1	62.9 - 83.4

Please refer to the package insert for the handling of XN CHECK BF.

Distributed by:
Sysmex America, Inc.
Lincolnshire, IL

350602
Rev. 3 25-Jan-23
May-24
2



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2024-08-11	Quality Control Data Due Date #1:	01-Jul-24	Quality Control Data Due Date #2:	11-Aug-24
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	41411301	41411302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.076	0.044 - 0.109	0.320	0.231 - 0.410
RBC-BF (10 ⁶ /μL)	0.025	0.020 - 0.030	0.078	0.067 - 0.089
TC-BF# (10 ³ /μL)	0.076	0.044 - 0.109	0.320	0.231 - 0.410
MN# (10 ³ /μL)	0.020	0.007 - 0.033	0.084	0.046 - 0.122
PMN# (10 ³ /μL)	0.056	0.026 - 0.087	0.236	0.160 - 0.311
MN%	26.1	10.2 - 42.0	26.4	15.6 - 37.2
PMN%	73.9	58.4 - 89.4	73.6	63.3 - 83.9

Please refer to the package insert for the handling of XN CHECK BF.