



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-10, XN-20) ANALYZERS

Expiration Date:	2024-09-08	Quality Control Data Due Date #1:	29-Jul-24	Quality Control Data Due Date #2:	08-Sep-24
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	41691301	41691302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.072	0.042 - 0.102	0.295	0.212 - 0.378
RBC-BF (10 ⁶ /μL)	0.024	0.019 - 0.028	0.073	0.063 - 0.083
TC-BF# (10 ³ /μL)	0.072	0.042 - 0.102	0.295	0.212 - 0.378
MN# (10 ³ /μL)	0.019	0.006 - 0.032	0.079	0.043 - 0.114
PMN# (10 ³ /μL)	0.053	0.024 - 0.081	0.216	0.147 - 0.286
MN%	26.8	10.4 - 43.1	26.7	15.7 - 37.6
PMN%	73.2	57.9 - 88.6	73.3	63.1 - 83.6
WBC ¹ (10 ³ /μL)	0.072	0.042 - 0.102	0.295	0.212 - 0.378

Please refer to the package insert for the handling of XN CHECK BF.

¹WBC for analysis in Blood Bank mode.

Distributed by:
Sysmex America, Inc.
Lincolnshire, IL

350602
Rev. 3 25-Jan-23
Jun-24
1



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-L ANALYZERS

Expiration Date:	2024-09-08	Quality Control Data Due Date #1:	29-Jul-24	Quality Control Data Due Date #2:	08-Sep-24
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	41691301	41691302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.073	0.043 - 0.104	0.298	0.215 - 0.381
RBC-BF (10 ⁶ /μL)	0.024	0.019 - 0.028	0.074	0.063 - 0.084
TC-BF# (10 ³ /μL)	0.073	0.043 - 0.104	0.298	0.215 - 0.381
MN# (10 ³ /μL)	0.020	0.008 - 0.032	0.081	0.044 - 0.117
PMN# (10 ³ /μL)	0.053	0.025 - 0.082	0.217	0.148 - 0.287
MN%	27.3	14.2 - 40.4	27.1	16.0 - 38.3
PMN%	72.7	57.4 - 88.0	72.9	62.7 - 83.1

Please refer to the package insert for the handling of XN CHECK BF.

Distributed by:
Sysmex America, Inc.
Lincolnshire, IL

350602
Rev. 3 25-Jan-23
Jun-24
2



XN CHECK™ BF

BODY FLUID CONTROL FOR SYSMEX XN-SERIES (XN-11, XN-21) ANALYZERS

Expiration Date:	2024-09-08	Quality Control Data Due Date #1:	29-Jul-24	Quality Control Data Due Date #2:	08-Sep-24
------------------	-------------------	--------------------------------------	-----------	--------------------------------------	-----------

Lot Number:	41691301	41691302
-------------	----------	----------

Control:	L1:Level 1	L2:Level 2
----------	------------	------------

PARAMETERS	MEAN	EXPECTED RANGE	MEAN	EXPECTED RANGE
WBC-BF (10 ³ /μL)	0.074	0.043 - 0.105	0.299	0.215 - 0.383
RBC-BF (10 ⁶ /μL)	0.024	0.020 - 0.029	0.073	0.063 - 0.084
TC-BF# (10 ³ /μL)	0.074	0.043 - 0.105	0.299	0.215 - 0.383
MN# (10 ³ /μL)	0.020	0.007 - 0.034	0.080	0.044 - 0.116
PMN# (10 ³ /μL)	0.054	0.025 - 0.083	0.219	0.149 - 0.289
MN%	27.3	10.7 - 44.0	26.8	15.8 - 37.8
PMN%	72.7	57.4 - 88.0	73.2	63.0 - 83.5

Please refer to the package insert for the handling of XN CHECK BF.